

QHerit™ Expanded Carrier Screen

Test Code: 94372 (X)

Specimen Requirements: Preferred: 6 mL (4 mL minimum) room-temperature whole blood; 1.5 mL (1 mL minimum) in each of 4 lavender-top (EDTA) or yellow-top (ACD) tubes

Alternatives: 10 to 15 mL room-temperature amniotic fluid (sterile plastic container) or 10 to 20 mg chorionic villi in 2 to 3 mL sterile saline or tissue culture medium (sterile plastic container)

Note: The QHerit™ Expanded Carrier Screen may not be ordered on fetal specimens. Fetal testing is available for each of the 22 individual diseases (ie, not as a panel).

CPT Codes*: 81200, 81205, 81209, 81220, 81242, 81243, 81250, 81251, 81255, 81257, 81260, 81290, 81330, 81400 (x4), 81401 (x3), 81479

CLINICAL USE

- Determine carrier risk for 22 clinically actionable, heritable diseases

CLINICAL BACKGROUND

Mendelian disorders occur before 25 years of age in an estimated 0.4% of live births.¹ Preconception and prenatal screening can identify carriers of genetic variants that cause these disorders and, thus, help couples plan and manage their pregnancies.

Carrier screening has traditionally targeted ethnicities that are at higher risk of disorders. An ethnicity-based approach presents difficulties for individuals who are multi-racial, adopted, or have incomplete or incorrect information about their ethnic backgrounds. This difficulty, combined with the advent of next-generation sequencing (NGS), has led to consideration of screening individuals regardless of ethnicity.

NGS has improved sequencing costs and efficiency, enabling the creation of expanded carrier screening (ECS) panels. These panels allow testing of many disorders at once and may improve outcomes and cost-effectiveness compared to traditional screening.² Studies have also shown that ECS panels are clinically valid for carrier screening^{3,4} and have clinical utility for infertile or at-risk couples.^{5,6}

The advantages of NGS have led to some published ECS panels that test for hundreds of conditions. This approach can lead to more than half of tested patients being identified as carriers,⁷ which can, in turn, lead to testing of partners and patient anxiety. Furthermore, testing for some conditions may not provide clinical utility (eg, conditions that have weak association with tested variants or are so rare that residual risk after a negative result is impossible to calculate).⁸ Thus, the conditions included on ECS panels should be carefully selected.

In March 2017, the American College of Obstetricians and Gynecologists (ACOG) released criteria for conditions included in ECS panels.⁷ The guidelines suggest that conditions meet “several” of the following criteria: 1) carrier frequency ≥ 1 in 100; 2) well-defined phenotype; 3) negative effect on quality of life; 4) cause of cognitive or physical disability; or 5) early onset. Furthermore, diagnosis before birth and disease treatment, management, or education should be possible.

In 2015, a joint statement about points to consider for ECS was also released by the following 5 medical organizations⁸: 1) the American College of Medical Genetics and Genomics (ACMG), 2) ACOG, 3) the National Society of Genetic Counselors (NSGC), 4) the Perinatal Quality Foundation, and 5) the Society for Maternal-Fetal Medicine (SMFM). The joint statement contains criteria that align with the 2017 ACOG guidelines. It also includes additional exclusion criteria. Specifically, conditions should not be included in ECS panels if: 1) testing cannot distinguish between adult and childhood onset; 2) associated alleles are relatively common in the general population but penetrance is low; or 3) other methods are more appropriate than molecular testing.

The QHerit Expanded Carrier Screen is consistent with the guidelines discussed above.^{7,8} The assay includes testing for clinically actionable variants in 24 genes related to 22 heritable diseases (**Table**) that are included in genetic testing guidelines from ACOG, ACMG, NSGC, or the Jewish Genetic Disease Consortium (JGDC).

INDIVIDUALS SUITABLE FOR TESTING

- Individuals of any ethnicity or geographic origin who are pregnant or considering pregnancy

Test Summary

Table. QHerit Expanded Carrier Screen: Tested Variants, Individual Tests, Related Guidelines, Detection Rates, and Residual Risk

Disease	Tested Gene(s)	Tested Variant(s)
Alpha-Thalassemia	<i>HBA1, HBA2</i>	-alpha3.7, -alpha4.2, -alpha20.5, --SEA, --MED, --FIL, --THAI, Constant Spring (c.427T>C)
Beta-Hemoglobinopathies (Including Sickle Cell Disease)	<i>HBB</i>	See footnote b.
Bloom Syndrome	<i>BLM</i>	2281del6/ins7 (c.2207_2212del(ATCTGA)insTAGATTC)
Canavan Disease	<i>ASPA</i>	IVS2-2A>G (c.433-2A>G), Y231X (c.693C>A), E285A (c.854A>C), A305E (c.914C>A)
Cystic Fibrosis	<i>CFTR</i>	See footnote d.
Dihydroipoamide Dehydrogenase Deficiency	<i>DLD</i>	Y35* (c.104dupA), G229C (c.685G>T)
Familial Dysautonomia	<i>IKBKAP</i>	R696P (c.2087G>C), IVS20+6T>C (c.2204+6T>C)
Familial Hyperinsulinism	<i>ABCC8</i>	IVS32-9G>A (c.3989-9G>A), F1387del (c.4160_4162delTCT)
Fanconi Anemia Type C	<i>FANCC</i>	IVS4+4A>T (c.456+4A>T), 322delG (c.67delG)
Fragile X Syndrome	<i>FMR1</i>	CGG triplet repeat number is reported
Gaucher Disease	<i>GBA</i>	IVS2+1G>A (c.115+1G>A), 84GG (c.84dupG), N370S (c.1226A>G), del55bp (c.1263_1317del55), V394L (c.1297G>T), D409H (c.1342G>C), L444P (c.1448T>C), R496H (c.1604G>A)
Glycogen Storage Disease Type IA	<i>G6PC</i>	R83C (c.247C>T), Q347X (c.1039C>T)
Joubert Syndrome 2	<i>TMEM216</i>	R73L (c.218G>T)
Maple Syrup Urine Disease	<i>BCKDHA</i>	R183P (c.548G>C), G278S (c.832G>A), E372X (c.1114G>T)
Mucopolidosis IV	<i>MCOLN1</i>	IVS3-2A>G (c.406-2A>G), 6.4kb_del (g.7586622_7593055del)
Nemaline Myopathy	<i>NEB</i>	2502bp del (c.7431+1917_7536+372del)
Niemann-Pick Disease Types A and B	<i>SMPD1</i>	L302P (c.911T>C), fsP330 (c.996delC), R496L (c.1493G>T), deltaR608 (c.1829_1831delGCC)
Spinal Muscular Atrophy	<i>SMN1, SMN2</i>	<i>SMN1</i> and <i>SMN2</i> copy numbers are reported
Tay-Sachs Disease	<i>HEXA</i>	Pseudodeficiency variants: R247W (c.739C>T), R249W (c.745C>T) Other variants: R178H/B1 variant (c.533G>A), G269S (c.805G>A), IVS9+1G>A (c.1073+1G>A), 1278TATC (c.1274_1277dupTATC), IVS12+1G>C (c.1421+1G>C), 7.6kb_del, Ex1 (c.2564_253+5128delinsG)
Usher Syndrome Type IF	<i>PCDH15</i>	R245* (c.733C>T)
Usher Syndrome Type IIIA	<i>CLRN1</i>	N48K (c.144T>G)
Walker-Warburg Syndrome	<i>FKTN</i>	F390fs (c.1167dupA)

ACMG indicates American College of Medical Genetics and Genomics; ACOG, American College of Obstetricians and Gynecologists; AJ, Ashkenazi Jewish; JGDC, Jewish Genetic Disease Consortium; NSGC, National Society of Genetic Counselors; SMFM, Society for Maternal-Fetal Medicine.

For footnotes, see Table Footnotes section in body text.

Individual Test (Test Code)	Guidelines That Indicate Testing	Ethnicity or Sex
Alpha-Globin Common Mutation Analysis (11175) ^a	ACOG (per ancestry) ^{7,9}	Mediterranean, Middle Eastern, Southeast Asian, African, Chinese, Asian Indian
Beta-Globin Complete (14974) ^c	ACOG (per ancestry) ^{7,9}	Mediterranean, Middle Eastern, Southeast Asian, African, Chinese, Asian Indian
Bloom Syndrome DNA Mutation Analysis (90872)	ACMG (AJ) ¹⁰ ; JGDC ¹¹	Ashkenazi Jewish
Canavan Disease Mutation Analysis (90905)	ACOG (AJ) ^{7,9} ; ACMG (AJ) ¹⁰ ; JGDC ¹¹	Ashkenazi Jewish Non-Ashkenazi Jewish
CFvantage [®] CF Expanded Screen (92068)	ACOG (all women) ^{7,9} ; ACMG (all women) ¹² ; NSGC (all women) ¹³ ; JGDC ¹¹	Ashkenazi Jewish Non-Hispanic Caucasian Hispanic American African American Asian American
Dihydroipoamide Dehydrogenase Deficiency (DLD Deficiency) (92046)	JGDC (lipoamide dehydrogenase deficiency [E3]) ¹¹	Ashkenazi Jewish
Familial Dysautonomia Mutation Analysis (90912)	ACOG (AJ) ⁹ ; ACMG (AJ) ¹⁰ ; JGDC ¹¹	Ashkenazi Jewish
Familial Hyperinsulinism (92045)	JGDC ¹¹	Ashkenazi Jewish
Fanconi Anemia DNA Mutation Analysis (90897)	ACMG (AJ) ¹⁰ ; JGDC ¹¹	Ashkenazi Jewish
XSense [®] , Fragile X w/Reflex (16313)	ACOG (family history) ^{7,9} ; ACMG (family history) ¹⁴ ; NSGC (family history) ¹⁵	Females
Gaucher Disease DNA Mutation Analysis (90907)	ACMG (AJ) ¹⁰ ; JGDC ¹¹	Ashkenazi Jewish
Glycogen Storage Disease Type IA Mutation Analysis (90915)	JGDC ¹¹	Ashkenazi Jewish Caucasian
Joubert Syndrome 2 (92050)	JGDC ¹¹	Ashkenazi Jewish
Maple Syrup Urine Disease (MSUD) Mutation Analysis (90909)	JGDC ¹¹	Ashkenazi Jewish
Mucopolipidosis IV Mutation Analysis (90899)	ACMG (AJ) ¹⁰ ; JGDC ¹¹	Ashkenazi Jewish
Nemaline Myopathy (92055)	JGDC ¹¹	Ashkenazi Jewish
Niemann-Pick Disease Mutation Analysis (90893)	ACMG (AJ) ¹⁰ ; JGDC ¹¹	Ashkenazi Jewish
SMA Carrier Screen (18041)	ACOG (family history) ^{7,9} ; ACMG (all women) ¹⁶ ; JGDC ¹¹	Caucasian Ashkenazi Jewish Asian African American Hispanic
Tay-Sachs Disease Mutation Analysis (90903)	ACOG (AJ, Cajun, French Canadian) ^{7,9} ; ACMG (AJ) ¹⁰ ; JGDC ¹¹	Ashkenazi Jewish French-Canadian General Population
Usher Syndrome Type IF (92047)	JGDC ¹¹	Ashkenazi Jewish
Usher Syndrome Type III (92048)	JGDC ¹¹	Ashkenazi Jewish
Walker-Warburg Syndrome (92051)	JGDC ¹¹	Ashkenazi Jewish

Detection Rate, %	Prior Carrier Risk	Residual Risk After a Negative Result		References for Rates and Risks
Up to 94	Varies by ethnicity	Reduced		17
99	Varies by ethnicity	Reduced		17 and footnote e
99	1/134	1/13,301		18
>97	1/55	<1/1,801		18
50	Not known			19,20
95	1/24	1/461		
90	1/25	1/241		
88	1/46	1/376		21-24
78	1/65	1/292		
53	1/94	1/199		
>95	1/107	<1/2,121		18,24
>99	1/31	<1/3,001		18
90	1/68	1/671		18
99	1/100	1/9,901		18
99	1/259	1/25,801		14,25,26
95	1/15	1/281		18
95	1/64	1/1,261		18
51	1/177	1/360		27
99	1/107	1/10,601		28
95	1/97	1/1,921		18
95	1/89	1/1,761		18
>95	1/168	<1/3,341		18
97	1/115	1/3,801		18
		2 SMN1 copy result	3 SMN1 copy result	
95	1/35	1/632	1/3,500	
90	1/41	1/350	1/4,000	29
93	1/53	1/628	1/5,000	
71	1/66	1/121	1/3,000	
91	1/117	1/1,061	1/11,000	
98	1/27	1/1,301		18
70	1/31	1/101		30
46	1/300	1/555		30
>75	1/147	<1/585		18
>95	1/120	<1/2,381		18
99	1/79	1/7,801		28

METHOD

- Fragile X syndrome (FXS): PCR and capillary electrophoresis; detect CGG repeat number
- Spinal muscular atrophy (SMA): real-time, allele-specific PCR; determine copy number of *SMN1* and *SMN2*
- Cystic fibrosis (CF): targeted multiplex NGS for 161 CF variants, including the 23 common variants recommended by ACOG/ACMG
- All other diseases: multiplex PCR of specific regions of 20 genes related to 19 diseases; NGS analysis follows PCR

INTERPRETIVE INFORMATION

Residual risk after a negative result for each gene/condition can be found in the **Table**. Results should be interpreted in conjunction with other laboratory and clinical findings.

FXS: Individuals with <45 CGG repeats are considered negative for FXS, while those with 45 to 54 CGG repeats are considered to be gray-zone allele carriers. No FXS-associated phenotype is expected for individuals with negative or gray-zone allele carrier results. If an expanded allele (>85 CGG repeats) is detected, both CGG repeat number and hypermethylation status are reported. Premutation allele carriers (repeat lengths of 55–200) are at risk of FXS-associated syndromes. Individuals with full mutation (>200 CGG repeats) and hypermethylation are predicted to be affected by FXS.

The associated risk of having a child with a premutation or a full mutation depends on the gender and mutation status of the parent and the gender of the child. This assay does not detect other mutations (eg, deletions and point mutations) that disrupt the function of the *FMR1* gene or protein.

SMA: A result of 1 or 0 *SMN1* copies indicates that the individual carries a disease-related deletion. This assay detects the copy number of *SMN1* and *SMN2*. This test cannot detect other pathogenic variants, nor can it identify silent carriers: individuals who have 2 or more copies of the

SMN1 gene on 1 chromosome and 0 copies of the *SMN1* gene on the opposite chromosome (eg, 2+0 carriers). The risk for pathogenic variants that cause SMA other than the deletions tested depends greatly on family history and clinical presentation. Therefore, if both members of a couple carry SMA deletions (with results indicating 1 or 0 copies of *SMN1* detected), or other pathogenic variants not detected by this assay, their children have a higher chance to be affected with SMA. This risk is 25% if both members of a couple each carry 1 *SMN1* deletion or other pathogenic variant.

All other genes/conditions: A positive result indicates that the individual carries a disease-related variant. Thus, the individual has an increased risk of having a child affected with the corresponding disease; testing of reproductive partners should be considered. If both members of a couple are carriers, their offspring have a 25% risk of being affected. The test report also indicates whether an individual has 1 copy of a pathogenic variant (carrier), 2 copies of the same variant (affected), or 1 copy each of 2 different variants (affected).

Additional assistance with interpretation of results is available for healthcare providers from our Genetic Counselors by calling Quest Genomics Client Services at 1.866.GENE.INFO (1.866.436.3463).

Visit QHerit.com/Clinicalinfo for additional clinical information and references.

References

1. Baird PA, Anderson TW, Newcombe HB, et al. Genetic disorders in children and young adults: a population study. *Am J Hum Genet.* 1988;42:677-693.
2. Azimi M, Schmaus K, Greger V, et al. Carrier screening by next-generation sequencing: health benefits and cost effectiveness. *Mol Genet Genomic Med.* 2016;4:292-302.
3. Lazarin GA, Haque IS, Nazareth S, et al. An empirical estimate of carrier frequencies for 400+ causal Mendelian variants: results from an ethnically diverse clinical sample of 23,453 individuals. *Genet Med.* 2013;15:178-186.

(continued on next page)

Table footnotes

^a QHerit Expanded Carrier Screen includes Hb Constant Spring, whereas the Alpha Globin Common Mutation Analysis test (test code 11175) does not.

^b *HBB* variants tested (including sickle cell disease and beta thalassemias): c.*111A>G, c.*110T>C, c.*96T>C, Hb D-Los Angeles (c.364G>C), Hb O-Arab (c.364G>A), c.321_322insG, c.316-2A>C, c.316-2A>G, c.316-3C>A, c.316-106C>G, c.316-125A>G, c.316-146T>G, c.316-197C>T, c.315+1G>A, c.287_288insA, c.251delG, c.230delC, c.216_217insA, c.203_204delTG, c.143_144insA, c.146_147insATCT, c.135delC, c.130G>T, c.126_129delCTTT, c.124_127delTTCT, c.118C>T, c.114G>A, c.112delT, c.93-1G>C, c.93-1G>A, c.93-21G>A, c.92+6T>C, c.92+5G>A, c.92+5G>C, c.92+5G>T, c.92+2T>A, c.92+2T>C, c.92+1G>A, c.92+1G>T, Hb Monroe (c.92G>C), c.92G>A, c.84_85insC, c.79G>T, HBE (c.79G>A), c.75T>A, c.59A>G, c.52A>T (LYS17*), c.51delC, c.48G>A, c.47G>A (Trp15), c.46delT, c.36delT, c.33C>A, c.27_28insG, c.25_26delAA, c.20delA, HBS (c.20A>T), HBC (c.19G>A), c.17_18delCT, c.2T>C, c.2T>G, c.1A>G, c.-78A>C, c.-78A>G, c.-79A>G, c.-80T>A, c.-81A>G, c.-136C>G, c.-137C>A, c.-137C>G, c.-137C>T, c.-138C>T, c.-138C>A, c.-140C>T, c.-151C>T.

^c QHerit Expanded Carrier Screen does not include full sequencing, whereas the Beta-Globin Complete test (test code 14974) does.

^d For *CFTR* variants tested, see online: QuestDiagnostics.com/testcenter/testguide.action?dc=TS_QHerit.

^e Based on Quest Diagnostics samples tested 2004-2016.

References (continued)

4. Bell CJ, Dinwiddie DL, Miller NA, et al. Carrier testing for severe childhood recessive diseases by next-generation sequencing. *Sci Transl Med*. 2011;3:65ra64.
5. Franasiak JM, Olcha M, Bergh PA, et al. Expanded carrier screening in an infertile population: how often is clinical decision making affected? *Genet Med*. 2016;18:1097-1101.
6. Martin J, Asan, Yi Y, et al. Comprehensive carrier genetic test using next-generation deoxyribonucleic acid sequencing in infertile couples wishing to conceive through assisted reproductive technology. *Fertil Steril*. 2015;104:1286-1293.
7. Committee on Genetics. Committee Opinion No. 690 Summary: Carrier Screening in the Age of Genomic Medicine. *Obstet Gynecol*. 2017;129:595-596.
8. Edwards JG, Feldman G, Goldberg J, et al. Expanded carrier screening in reproductive medicine—points to consider: a joint statement of the American College of Medical Genetics and Genomics, American College of Obstetricians and Gynecologists, National Society of Genetic Counselors, Perinatal Quality Foundation, and Society for Maternal-Fetal Medicine. *Obstet Gynecol*. 2015;125:653-662.
9. Committee on Genetics. Committee Opinion No. 691: Carrier Screening for Genetic Conditions. *Obstet Gynecol*. 2017;129:e41-e55.
10. Gross SJ, Pletcher BA, Monaghan KG. Carrier screening in individuals of Ashkenazi Jewish descent. *Genet Med*. 2008;10:54-56.
11. Jewish Genetic Diseases Consortium. Jewish Genetic Diseases. 2014; http://www.jewishgeneticdiseases.org/wp-content/uploads/2014/11/2014_2_JGDC_Diseases.pdf. Updated 2014. Accessed June 6, 2017.
12. Watson MS, Cutting GR, Desnick RJ, et al. Cystic fibrosis population carrier screening: 2004 revision of the American College of Medical Genetics mutation panel. *Genet Med*. 2004;6:387-391.
13. Langfelder-Schwind E, Karczeski B, Strecker MN, et al. Molecular testing for cystic fibrosis carrier status practice guidelines: recommendations of the National Society of Genetic Counselors. *J Genet Couns*. 2014;23:5-15.
14. Monaghan KG, Lyon E, Spector EB. ACMG Standards and Guidelines for fragile X testing: a revision to the disease-specific supplements to the Standards and Guidelines for Clinical Genetics Laboratories of the American College of Medical Genetics and Genomics. *Genet Med*. 2013;15:575-586.
15. Finucane B, Abrams L, Cronister A, et al. Genetic counseling and testing for FMR1 gene mutations: practice guidelines of the National Society of Genetic Counselors. *J Genet Couns*. 2012;21:752-760.
16. Prior TW, Nagan N, Sugarman EA, et al. Technical standards and guidelines for spinal muscular atrophy testing. *Genet Med*. 2011;13:686-694.
17. Origa R, Moi P. Alpha-Thalassemia. <https://www.ncbi.nlm.nih.gov/books/NBK1435/>. Updated December 29, 2016. Accessed June 6, 2017.
18. Scott SA, Edelmann L, Liu L, et al. Experience with carrier screening and prenatal diagnosis for 16 Ashkenazi Jewish genetic diseases. *Hum Mutat*. 2010;31:1240-1250.
19. Kaul R, Gao GP, Matalon R, et al. Identification and expression of eight novel mutations among non-Jewish patients with Canavan disease. *Am J Hum Genet*. 1996;59:95-102.
20. Sistermans EA, de Coo RF, van Beerendonk HM, et al. Mutation detection in the aspartoacylase gene in 17 patients with Canavan disease: four new mutations in the non-Jewish population. *Eur J Hum Genet*. 2000;8:557-560.
21. Friedman KJ, Leigh MW, Czarnecki P, et al. Cystic fibrosis transmembrane-conductance regulator mutations among African Americans. *Am J Hum Genet*. 1998;62:195-196.
22. Heim RA, Sugarman EA, Allitto BA. Improved detection of cystic fibrosis mutations in the heterogeneous U.S. population using an expanded, pan-ethnic mutation panel. *Genet Med*. 2001;3:168-176.
23. Schrijver I, Ramalingam S, Sankaran R, et al. Diagnostic testing by CFTR gene mutation analysis in a large group of Hispanics: novel mutations and assessment of a population-specific mutation spectrum. *J Mol Diagn*. 2005;7:289-299.
24. Sugarman EA, Rohlfes EM, Silverman LM, et al. CFTR mutation distribution among U.S. Hispanic and African American individuals: evaluation in cystic fibrosis patient and carrier screening populations. *Genet Med*. 2004;6:392-399.
25. Dombrowski C, Levesque S, Morel ML, et al. Premutation and intermediate-size FMR1 alleles in 10572 males from the general population: loss of an AGG interruption is a late event in the generation of fragile X syndrome alleles. *Hum Mol Genet*. 2002;11:371-378.
26. Rousseau F, Rouillard P, Morel ML, et al. Prevalence of carriers of premutation-size alleles of the FMRI gene—and implications for the population genetics of the fragile X syndrome. *Am J Hum Genet*. 1995;57:1006-1018.
27. Chou JY, Mansfield BC. Mutations in the glucose-6-phosphatase-alpha (G6PC) gene that cause type Ia glycogen storage disease. *Hum Mutat*. 2008;29:921-930.
28. Hoffman JD, Park JJ, Schreiber-Agus N, et al. The Ashkenazi Jewish carrier screening panel: evolution, status quo, and disparities. *Prenat Diagn*. 2014;34:1161-1167.
29. Hendrickson BC, Donohoe C, Akmaev VR, et al. Differences in SMN1 allele frequencies among ethnic groups within North America. *J Med Genet*. 2009;46:641-644.
30. Monaghan KG, Feldman GL, Palomaki GE, et al. Technical standards and guidelines for reproductive screening in the Ashkenazi Jewish population. *Genet Med*. 2008;10:57-72.

* The CPT codes provided are based on AMA guidelines and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.

QuestDiagnostics.com

Quest, Quest Diagnostics, QHerit, CFvantage, XSense, any associated logos, and all associated Quest Diagnostics registered or unregistered trademarks are the property of Quest Diagnostics. All third-party marks—® and ™—are the property of their respective owners. © 2017 Quest Diagnostics Incorporated. All rights reserved. TS6999 06/2017